

# Product parameters

Product name	Silicone Facial Mask with Cable Facial Beauty Therapy Device for Skin Rejuvenation   Dongguan OEM Manufacturer
Cable length	120cm or customized
Cable color	white black or customized
Function	TENS EMS Transmission/Massager ECG/EKG/EEG/EMG Adapter Cable
Minimum order quantity	100pcs

Manufactured by a professional ISO13485-certified OEM factory in Dongguan, China, this Silicone Facial Mask with Cable Beauty Therapy Device is engineered for high-performance facial care. Similar to our [tVNS Ear Clip Electrode for Auricular Vagus Nerve Stimulation](#), this device demonstrates our expertise in developing advanced medical-grade and beauty therapy equipment.

The mask features medical-grade silicone for comfort, durability, and biocompatibility, ensuring safe use for sensitive skin. With integrated conductive technology, the device transmits low-frequency microcurrents or light therapy through the cable connection, promoting collagen production, tightening facial muscles, and improving skin elasticity.

Ideal for beauty salons, spas, and home skin care routines, this device helps reduce fine lines, restore hydration, and enhance skin radiance. Its ergonomic design ensures full face coverage while maintaining flexibility for user comfort. For complementary equipment, explore our [Customized 10 Pin Pogo Pin Magnetic Connector to 10 Leads 4.0mm Female Snap ECG Cable](#), which highlights our precision engineering in cable manufacturing.

We also provide OEM/ODM services, enabling custom branding, packaging, and design to meet diverse market needs. With international certifications such as CE and RoHS, this facial mask device is a trusted choice for global distributors and skin care brands.

## Product photos:










## Company information. Overview



## About Us

10+ Years Experienced Professional Manufacturer  
of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables,  
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

# FAQ

## **Q1: Can I start a sample order before ordering more?**

A: Of course, Heymed provides free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

## **Q2: What payment methods do you accept?**

A: PayPal, T/T, bank transfer, Western Union and L/C

## **Q3: How does your factory do in terms of quality control?**

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from the beginning to the end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

## **Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?**

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

## **Q5: Delivery time and warranty? What shipping service will you use to ship my product?**

A: Our delivery time is about 2-4 weeks, and we can provide a 12-month warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

## **Q6: Where is your factory located? How can I visit there?**

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China