

Product Details:

The [Reliable Type C to Pulse Oximeter Clip Cable](#) is an essential component for healthcare professionals seeking efficient and stable data transmission between pulse oximeters and monitoring devices. This high-quality cable is specifically designed to ensure seamless connectivity in critical medical environments, such as hospitals, clinics, and emergency care settings.



Applications:

- **Hospitals and Clinics:** The cable is ideal for use in hospital settings, where rapid and accurate monitoring of patients' oxygen saturation levels is crucial. It connects pulse oximeters to monitoring systems, providing healthcare providers with real-time data.
- **Emergency Medical Services (EMS):** In emergency situations, where every second counts, this cable ensures fast and reliable data transmission from pulse oximeter devices to monitoring equipment, allowing paramedics to make quick, informed decisions.
- **Home Healthcare:** It is also suitable for home healthcare settings, where patients require

continuous monitoring of vital signs, including blood oxygen levels, and caregivers need a reliable cable for medical devices.

Key Features:

- **Premium Build Quality:** The cable features a high-quality [Type C](#) connector, offering fast and secure data transmission. The clip end is designed for a perfect fit with pulse oximeter sensors, ensuring a stable and precise connection.
- **Medical-Grade Materials:** Crafted from medical-grade **silicone** and **plastic**, this cable is non-toxic, hypoallergenic, and resistant to frequent cleaning and sterilization, which is essential in medical environments.
- **Durable and Long-Lasting:** Built to last with an estimated lifespan of **3,000+ insertions**, it is highly durable and resistant to wear and tear, minimizing the need for frequent replacements.
- **Flexible and Tangle-Free:** With its flexible design, this cable can be easily maneuvered and positioned as needed. Its tangle-free construction prevents frustration and saves time, which is crucial in fast-paced medical situations.
- **Compatible with Various Medical Devices:** While primarily designed for pulse oximeters, this cable is compatible with other devices like [ECG cables and electrode cables](#), making it a versatile addition to your medical equipment inventory.
- **Plug-and-Play Simplicity:** No complicated setup is required. Simply plug the **Type C** connector into your device and the clip end into the pulse oximeter sensor. The cable is easy to use, allowing healthcare professionals to focus on patient care without technical distractions.



Certifications and Compliance:

This cable meets rigorous standards to ensure the highest levels of quality and safety in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Guarantees the product is manufactured under a robust quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration, meeting their standards for medical devices.
- **CE MDR Certified:** Complies with the European Union's Medical Device Regulation (MDR), ensuring safety and effectiveness for medical use.


Why Choose Our Type C to Pulse Oximeter Clip Cable?

- **Reliability in Critical Care:** This cable ensures continuous, accurate data transmission, providing healthcare providers with real-time, reliable readings of blood oxygen saturation levels.
- **Versatility:** With broad compatibility across various medical devices, including **ECG monitors** and

electrode cables, it offers a unified solution for different types of monitoring equipment.

- **Cost-Effective:** The high durability and long lifespan of the cable reduce the frequency of replacements, offering significant cost savings over time.
- **Improved Efficiency:** The cable's flexibility, durability, and ease of use allow healthcare professionals to focus on delivering the best patient care rather than dealing with equipment issues.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed provides free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from the beginning to the end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide a 12-month warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China