


|                    |                                                                                                                |
|--------------------|----------------------------------------------------------------------------------------------------------------|
| Product name       | USB Type C to 4.0mm Female Snap Cable with 2 Leads Wire   Custom Medical Lead Wire from Dongguan Cable Factory |
| application        | Stimulation machine, physical therapy instrument                                                               |
| length             | 10cm, 20cm, 50cm, 100cm or customized                                                                          |
| color              | Gray, white, black and customized                                                                              |
| Conductor material | Tinned copper                                                                                                  |
| Package            | Individually packaged                                                                                          |

This Type C to 4.0mm female snap cable is an innovative solution for connecting modern USB-C enabled medical or therapy devices with traditional snap electrode pads. Manufactured in Dongguan, China, this cable is ideal for TENS, EMS, and other biofeedback systems, offering custom options and reliable quality from a trusted Chinese factory and supplier.

Manufactured by a professional factory based in Dongguan, China, we are an ISO13485 and ISO9001 certified supplier with rich experience in OEM and ODM production. This product is widely used in TENS therapy, EMS muscle stimulators, and other electrotherapy equipment. For similar applications, we also offer a [USB-C to 4.0mm female snap electrode cable with silicone earhook](#), which ensures reliable signal transmission and enhanced patient comfort, as well as a [custom 6.35mm audio plug to dual 2.0mm electrode pin cable](#) designed for professional physiotherapy settings.

Custom lengths, packaging, and private labeling are available upon request.





### **Certification and Compliance:**

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

### **Company information. Overview**



## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10 years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



# FAQ

## **Q1: Can I start a sample order before ordering more?**

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

## **Q2: What payment methods do you accept?**

A: PayPal, T/T, bank transfer, Western Union and L/C

## **Q3: How does your factory do in terms of quality control?**

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

## **Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?**

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

## **Q5: Delivery time and warranty? What shipping service will you use to ship my product?**

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

## **Q6: Where is your factory located? How can I visit there?**

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China